

IV-STY[®]

" Innovative design Clip that secures IV drip lines whilst removing the risk of needle-stick injuries"



Sam McDonald - Nurse Practitioner and Inventor of IV-Sty[®]

The IV-Sty[®] - is a simple disposable clip that securely connects intra-venous lines enabling two sources of fluids to flow into one tube

Easier for the practitioner – Clips onto many IV line systems, removes the improvised use of clinical tape, improves line of sight to the "Piggy Back" port to ensure the free flow of fluid to patients.

Hygienic – Simply disposed of after use, removes solutions that can create infection control issues such as leaving residue adhesive to the line.

Health and Safety – Prevents the accidental disconnection of lines or needle becoming detached, thus reducing the probability of needle-stick injuries to staff or patients.

Cost effective – Removes the likelihood of damage/waste to secondary lines when changing.

The IV-Sty[®] has been designed to improve patient and practitioner safety to ensure 'best practice' when delivering secondary IV fluids.



Licensed to MDTi by Cornwall & Isles of Scilly NHS Primary Care Trust

EASY ~ EFFICIENT ~ EFFECTIVE

MDTi, The Kace Building, Victoria Passage, Wolverhampton WV1 4LG
Tel: +44 (0) 1902 778380 Fax: +44 (0) 1902 421360
Web: www.mdti.co.uk Email: info@mdti.co.uk

IV-Sty[®]
an

MDTi
product

OTHER PRODUCTS ASSOCIATED WITH THE IV-STY[®]



THE HOO-KON[®] PORTAL IV STAND

The Hoo-kOn[®] meets health and safety and best infection control practice, approved fit-for-purpose and developed within the NHS to deliver improved patient care in the management of IV/fluid therapy.

Benefits:-

- Cost effective
- All parts are easy to clean and are autoclavable
- Meets infection control requirements
- Easy to store away and can be carried in any small vehicle
- Lightweight and durable



NHS Supplychain Code: FAM090

The Hoo-kOn[®] has been designed for Community and Emergency use when dealing with patients that require IV/fluid therapy. The Hoo-kOn[™] comes in a modular form for ease of storage and cleaning whilst also lightweight so as to be carried from one site to another easily.



HOOK-IT[®] RANGE OF IV HOOKS

The HOOK-IT[®] family of drip hooks are simple to use products replacing many of the uses of traditional IV drip hooks and offer the following benefits:-

- Increases space round the patient's bed (no need for catheter stands)
- No risk of tripping or falling over catheter stands
- Easily cleaned with soap and water or alcohol and/or can be steamed sterilised.
- Can easily be moved from one area to another without the risk of transferring infection.
- Cost only a fraction of the price of catheter stands.



'C HOOKS' - CATHETER

NHS Supplychain Code: FKW006

Important: things to know about the Hook-It[®]

The HOOK-IT[®] range of hooks is manufactured from high quality materials that withstand steam sterilisation. Hooks can withstand multiple sterilisation cycles but MDTi recommends that they are replaced after two years of use.

Important notes on the Hoo-kOn[®] Portal IV Stand

The Hoo-kOn[®] stand has been designed to allow for ease of assembly in less than 2 minutes and is able to be assembled in any required location quickly. Its lightweight design allows for it to be picked up and placed in a position that allows IV/Fluid therapy to be given effectively. The stand can be mounted with an IV Pump, the Hoo-kOn[®] stand will hold an IV pump up to a weight of 6.7 kgs along with additional 2ltr fluid bag.

~ SIMPLE COST EFFECTIVE SOLUTIONS ~
CLINICAL INNOVATIONS.

Enhancing safety and efficiency

Helping to prevent the spread of MRSA

MDTi, The Kace Building, Victoria Passage, Wolverhampton WV1 4LG

Tel: +44 (0) 1902 778380 Fax: +44 (0) 1902 421360

Web: www.mdti.co.uk Email: info@mdti.co.uk